

Recommended protocol for ECAD

Obtained in run 53

07 Mar 2018

Immunostainer

Type: Ventana Benchmark XT

Primary antibody

Clone: GM016
Producer: Genemed
Product no. / lot no.: 61-0028 / 30119548
Diluent: Emerald - Antibody Diluent
Diluent producer / no.: Cell Marque / 936B
Dilution factor: 1:150
Incubation time / temperature: 44 min. / 37°C

Epitope retrieval, HIER

Device: On Board / On Machine
Buffer: Ventana CC1
Heating time at max. temp.: 60 min.
Maximum heating temp.: 100°C

Visualization system

Producer: Ventana
Product / no: UltraView Universal DAB Detection Kit / 760-500
Incubation time polymer: 12 min.
Incubation temperature: 37°C

Chromogen

Producer: Ventana
Product / no: ultraView Universal DAB Detection Kit / 760-500
Incubation time / temperature: 8 min. / 37°C
Enhancement: CuSO₄

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.