

## Recommended protocol for SOX10

Obtained in run 60

13 Aug 2020

### Immunostainer

Type: Dako Autostainer

### Primary antibody

Clone: ZM10  
Producer: Zeta Corportion  
Product no. / lot no.: Z2293 / 1511062F  
Diluent: Normal Antibody Diluent  
Diluent producer / no.: Zeta Corporation / ZD19  
Dilution factor: 1:200  
Incubation time / temperature: 45 min. / 24°C

### Epitope retrieval, HIER

Device: Pressure Cooker  
Buffer: Zeta Corporation Tris-EDTA HIER Solution pH 9.0  
Heating time at max. temp.: 45 min.  
Maximum heating temp.: 95°C

### Visualization system

Producer: Zeta Corporation  
Product / no: Zeta Universal HRP Polymer / ZD10  
Linker: None  
Incubation time polymer: 30 min.  
Incubation temperature: 24°C

### Chromogen

Producer: Zeta Corporation  
Product / no: DAB Chromogen / ZD11  
Incubation time / temperature: 5 min. / 24°C  
Enhancement: None

### Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.