

Recommended protocol for MLH1

Obtained in run 49

28 Dec 2016

Immunostainer

Type: LabVision 480

Primary antibody

Clone: ES05
Producer: Dako
Product no. / lot no.: M3640 / 10085035
Diluent: Antibody Diluent
Dilution factor: 1:50
Incubation time / temperature: 60 min. / 24°C

Epitope retrieval, HIER

Device: Microwave oven
Buffer: Dako Target Retrieval solution
Heating time at max. temp.: 24 min.
Maximum heating temp.: 99°C

Visualization system

Producer: IL Immunologic
Product / no: BrightVision / DPVO999HRP
Linker: None
Incubation time polymer: 30 min.
Incubation temperature: 24°C

Chromogen

Producer: Dako
Product / no: DAB+ / K5007 / K4065
Incubation time / temperature: 10 min. / 24°C
Enhancement: CuSO₄

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.