

Recommended protocol for PMS2

Obtained in run 53

26 Feb 2018

Immunostainer

Type: Thermo Autostainer 36/48/72

Primary antibody

Clone: A16-4
Producer: BD Biosciences
Product no. / lot no.: 556415 / 7103847
Diluent: Normal Antibody Diluent
Dilution factor: 1:100
Incubation time / temperature: 60 min. / 20°C

Epitope retrieval, HIER

Device: PT-link / PT-module
Buffer: Std Tris-EDTA / EGTA pH 9
Heating time at max. temp.: 10 min.
Maximum heating temp.: 96°C

Visualization system

Producer: Immunologic
Product / no: Brightvision Poly-HRP-Anti MS/Rb/Rt IgG / VWRKDPVO9
Linker: None
Incubation time polymer: 30 min.
Incubation temperature: 20°C

Chromogen

Producer: Immunologic
Product / no: Bright-DAB Substrate / VWRKBS04-999
Incubation time / temperature: 7 min. / 20°C
Enhancement: None

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.