

Recommended protocol for BSAP

Obtained in run

23 Sep 2015

Immunostainer

Type: Manual

Primary antibody

Clone: DAK-Pax5
Producer: Dako
Product no. / lot no.: M7307 / 00049202
Diluent: Antibody Diluent
Dilution factor: 1:50
Incubation time / temperature: 30 min. / 20°C

Epitope retrieval, HIER

Device: PT-link / PT-module
Buffer: Dako TRS High pH (3-1)
Heating time at max. temp.: 30 min.
Maximum heating temp.: 96°C

Visualization system

Producer: Dako
Product / no: Real EnVision Rabbit/Mouse / K5007
Incubation time polymer: 30 min.
Incubation temperature: 20°C

Chromogen

Producer: Dako
Product / no: DAB+ / K5007
Incubation time / temperature: 5 min. / 20°C
Enhancement: None

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.