

## Recommended protocol for MLH1

Obtained in run

26 Sep 2015

### Immunostainer

Type: Autostainer 480, LabVision

### Primary antibody

Clone: ES05  
Producer: Leica/novocastra  
Product no. / lot no.: NCL-L-MLH1 / 6013079  
Diluent: NS  
Dilution factor: 1:100  
Incubation time / temperature: 60 min. / 20°C

### Epitope retrieval, HIER

Device: Microwave oven  
Buffer: Std Tris-EDTA / EGTA pH 9  
Heating time at max. temp.: 15 min.  
Maximum heating temp.: 100°C

### Visualization system

Producer: LabVision  
Product / no: NS / TL-125-QHL  
Linker: NS  
Incubation time linker: 10 min.  
Incubation time polymer: 10 min.  
Incubation temperature: 20°C

### Chromogen

Producer: LabVision  
Product / no: DAB / TA-125-QHSX  
Incubation time / temperature: 10 min. / 20°C  
Enhancement: None

### Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.